

Evaluating biomass reactions with cosolvents/CO₂ under subcritical and supercritical conditions

Sipho C. Ndlela, Norman K. Olson

Iowa Energy Center, Iowa State University, BECON Facility, 1521 West F Avenue, Nevada, IA 50201, USA

INTRODUCTION

Cellulose and corn stover were subjected to a reaction with polar cosolvents (methanol and water)/CO₂ at subcritical and supercritical conditions. Resulting products from a cellulose reaction at 300°C analyzed by GCMS consisted of a mixture of organic acids, sugars, and several other components. Examples of organic acids included acetic, hexanoic, heptanoic and malonic acids, whereas sugars were levoglucosan, D-Allose, and D-Galactose. Furthermore product samples from a corn stover reaction at 300°C suggested presence of ketones, aromatic hydrocarbons, and methyl esters. As evidenced by these results, reactions at high pressure are plausible although they are accompanied by a complex product composition, and an unreacted residue of cellulose and lignin.

MATERIALS & METHODS

Cellulose (20 micron) characterization: (TGA, SEM, XRD) TGA: (a.) Non-isothermal conditions: 50 to 900°C, 5°C/min, N₂

- (b.) Isothermal conditions: 200°C to 400°C, 200min., N₂ SEM: Hitachi S-2460 N variable pressure, and SEM X-ray
- spectrometer XRD: Siemens D 500 x-ray diffractometer with CuK,



Cellulose Characterization Results

In Figure 1 weight loss at 85°C associated with loss of physisorbed water, at 305°C associate with loss of O. There were no significant SEM morphological changes under TGA conditions.

- In Figure 2, no weight loss at 200°C (grey/white color), gradual weight loss at 250°C (grey color), accelerated weight loss at 300 and 400 °C (black/ grey color)
- SEM x-ray spectrometry data:
- TGA at 200°C looks like the reference cellulose material (Fig 4a).

TGA at 300°C, 400°C and -900°C are similar and have lost O compared to the reference sample (Fig 4b).

MATERIALS & METHODS

Cellulose and corn stover experiments using a high pressure batch and flow reactor. Cosolvents: CO₂, Water, Methanol, Analysis instruments: SEM: Hitachi S-2460 N variable Pressure and GCMS, Agilent GC model 6890 integrated with a micromass GCT, DB5 column, Ei and Ci (NH₃ and CH₄)



 Fig.6a Batch reactor product
 Fig.6b Flow reactor

 (L-R): Rotavaped, Crude and
 product crude (L-R): at

 recovered methanol (300°C)
 250 °C and 300 °C



Figure 7. GC trace from analysis of RTVPcrude sample after a reaction with methanol at 300°C and 2500 psi (Top: cellulose, Bottom: corn stover)

Cellulose-Corn Stover cosolvent reaction results summary

- GCMS is a useful analytical tool for evaluating the product sample post high pressure cellulose/corn stover-methanol reactions.
- CO₂ addition during a high pressure reaction had a minimal effect in converting cellulose in the temperature range 200 to 300°C.
- Decreased amounts of cellulose residue remained after reactions above 300°C.
- Samples from batch reactor were much darker due to extended heating and cooling period when compared to corresponding conditions with a continuous flow reactor.
- Resulting products from a cellulose reaction at 300°C analyzed by GCMS consisted of a mixture of organic acids, sugars, and several other components. Examples of organic acids included acetic, hexanoic, heptanoic and malonic acids, whereas sugars were levoglucosan, D-Allose, and D-Galactose. Furthermore product samples from a corn stover reaction at 300°C suggested presence of ketones, aromatic hydrocarbons, and methyl esters. As evidenced by these results, reactions at high pressure are plausible although they are accompanied by a complex product composition, and an unreacted residue of cellulose and lignin.

ACKNOWLEDGMENTS

lowa Energy Center, DOE Fund: DE-FG36-06GO86014

Dr. Brent Shanks, (Iowa State University Chemical and Biological Engineering Department)

Dr. Ronald Holtan, (Iowa State University Chemistry Department)

Dr. Walter Trahanovsky, (Iowa State University Chemistry Department)